



**Medrio R42.9 Minor Release Notes
2025-07-12**

Table of Contents

Introduction (3)

- Release Notes Overview (3)
- Release Preparation (3)
- What's New (3)

Medrio EDC (4)

- Monitoring UI Update (4)
- Log In - Provide the Ability to Display a Customer Agreement upon Log In to Medrio (6)
- Add Study Title to User Audits (8)
- Ability to Download Previously Attached Files in Audit Trail (9)

Medrio ePRO (10)

- Localization - Danish Language (10)

Introduction

Release Notes Overview

Medrio's R42.9 release is nearly here! The release notes below are intended to provide you with specific details on each of the updates, including an in-depth feature description, why the update was implemented, how to use the update, and the impact on existing and future studies.

Release Preparation

There is no webinar for this release. For all information related to R42.9, visit the Medrio Community at: [Release Notes](#).

If you have any questions or require assistance accessing the Medrio Community, please submit a support ticket by clicking **Help > Contact Us > Submit a Ticket** in Medrio.

What's New:


The R42.9 minor release includes updates to Medrio **CDMS/EDC** and **eCOA/ePRO**. CDMS/EDC updates strengthen study build efficiency, improve audit readiness, and enhance data traceability. eCOA/ePRO updates improve participant experiences.

Continue reading below for specific details regarding R42.9 changes.

Note: *Fixed and Resolved Issues will be posted on 2025-07-14 to the Medrio Community.*

Medrio EDC Monitoring UI Update

Description: The monitoring configuration page has been updated with a new user interface for selecting the forms and variables that are required for monitoring.




Context: This update is in preparation for future monitoring enhancements.

Activation: No action needed - this feature is turned on automatically for all new and existing studies.

Impact		
Product	Availability	Risk
EDC: Monitoring UI	Available for existing studies and studies provisioned after release	Low - this is simply a UI update to prepare for future enhancements to monitoring.

Log In - Provide the Ability to Display a Customer Agreement upon Log In to Medrio

Description: When Medrio updates our End-User License Agreement (EULA), users will now be prompted to agree to the latest version of the EULA before they can continue using Medrio. To make viewing the agreement easier, we have also added a User Agreement option to the Help menu.



Context: Previously, Medrio had no way of updating our EULA without a code change.

Activation: No action needed - this is automatically enabled for new and existing studies.

Impact		
Product	Availability	Risk
EDC: EULA	Available for existing studies and studies provisioned after release	Low – this simply asks the customer to agree to terms for using Medrio.

Add Study Title to User Audits

Description: The User Audit Log that can be downloaded via User Administration in the Manage tab now includes the Study Title in the file name and in a new Study Title column in the audit log.

Context: The study title was previously missing from the User Audit Log, making it difficult to pinpoint which study the User Audit Log was downloaded from.

Activation: No action needed - this is automatically enabled for new and existing studies.

Impact		
Product	Availability	Risk
EDC: User Audit Log	Available for existing studies and studies provisioned after release	Low – this feature simply adds the study title to the data and the file name for the User Audit Log. It does not touch the actual audit data.

Ability to Download Previously Attached Files in Audit Trail

Description: Files that were previously attached to file type variables and were subsequently replaced by other files can now be downloaded from the form-level audit logs. The Audit entry for each file now includes the file name, which is linked to the download of the file. You now have a more comprehensive document history for any file uploads that are versioned or replaced over the course of a study.

Context: Prior to this feature, previously attached files could not be downloaded at all, preventing auditors from being able to see the history of the attached files. This feature increases the accessibility of historical data by making previously attached and then replaced files downloadable. This does not alter existing data but provides a more complete audit trail and access to prior versions of uploaded files.

Activation: No action needed - this is automatically enabled for new and existing studies.

Impact		
Product	Availability	Risk
EDC: Audit Log	Available for existing studies and studies provisioned after release	Low – this feature provides access to previously existing files that were already stored. It's not modifying current data or how new data is captured. It's a read/download operation on historical versions.

Medrio ePRO

Localization - Danish Language

Description: Danish is now available for language selection in the participant portal and in the Add New Subject Information pop-up box.

The elements that display in Danish include:

- Navigation buttons
- Web page text
- Action buttons and other translated non-configurable UI elements (both on the remote and in-clinic pages)

Context: Participants were unable to complete ePROs in the Danish language.

Activation: No action needed - the feature is automatically available for studies using ePRO.

Impact		
Product	Availability	Risk
ePRO: Localization	Available for existing studies and studies provisioned after release	Low – this adds the ability to display system labels, buttons, and messages in the local language.